

# Anti-Sf9 HCP Antibody

Product No.: 1301312-Ab01

## ■Product information

<b>Size</b>	1 mg
<b>Species Reactivity</b>	<i>Spodoptera frugiperda</i>
<b>Host/Isotype</b>	Sheep / IgG
<b>Class</b>	Polyclonal
<b>Type</b>	Antibody
<b>Conjugate</b>	Unconjugated
<b>Immunogen</b>	Sf9 HCPs produced by sonication cell lysis
<b>Purification</b>	Affinity Purified
<b>Buffer</b>	PBS
<b>Contains</b>	No preservatives
<b>Applications</b>	HCP Coverage evaluations, Western Blot, ELISA
<b>Storage</b>	Maintain at 2-8°C for short-term use. For long-term storage, please store below -65°C.
<b>Validity period</b>	12 months

## ■Attentions

This product is for research use only and is not intended for diagnostic purposes.

## ■Related Product

Product Name	Product No.	Quantity
Sf9 HCP ELISA Kit (One-step ELISA)	1301312	96 tests/Kit

Effective date: 22 Jul. 2025

## Support & Contact

**SHENTEK**

Huzhou Shenke Biotechnology Co., Ltd.

[www.shentekbio.com](http://www.shentekbio.com)

Address: 8th Floor, 6B Building, No.1366 Hongfeng Road, Huzhou313000, Zhejiang Province, China

E-mail: [info@shentekbio.com](mailto:info@shentekbio.com)

Phone: +1 (908) 822-3199 / (+86) 400-878-2189